



Eyepieces and Monitors

What needs to be cleaned?

- ❑ Eyepieces on microscopes (6 in Plant Lab, 1 in Cloud Forest)
- ❑ Eyepieces on refractometers (2 in Plant Lab)
- ❑ Eyepieces on kaleidoscopes (6 along handrail in Rain Forest)

How do we clean them?

1. Pick up the small white bucket of eyepiece cleaning supplies at the Information Desk. It should contain: **q-tips, a glass cloth, a spray bottle of lens cleaner, a spray bottle of screen cleaner, and the key to the microscope cases.**
2. Spray lens cleaner onto a clean q-tip. (Don't spray the eyepiece of the equipment directly, since it might flood the lens).
3. Wipe the damp q-tip on the eyepiece, cleaning the area around the lens where people's eyes might touch it. The goal is to clean off germs and remove eye makeup.
4. To clean monitors, spray screen directly with **screen cleaner**. Use glass cloth to wipe down.